

Anti-gp120(92RW020)(Clade A)



Product Name	Anti- gp120(92RW020)(Clade A)
Catalog No	IT-001-001M1
Source	Mouse monoclonal antibody, clone 12G12
Isotype	IgG2a
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing gp140(92RW020)(Clade A) protein (a.a. 24-680) (Genebank accession AAT67478.1).
Applications	ELISA, WB, IP, IF etc.
Epitope	Linear
Specificity	Reacts with HIV-1 gp120(92RW020)(Clade A) or related strains only as tested so far. Minimal cross-reactivity to other subtypes.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

ELISA Profile of IT-001-001M1

