

Anti-gH/gL Complex(EBV/B95-8)



Product Name	Anti- gH/gL Complex(EBV/B95-8)
Catalog No	IT-005-030M9
Source	Mouse monoclonal antibody, clone 7E7
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	Co-immunization with DNA vaccine expressing the ectodomains of gH protein (a.a.1-529) (Genebank #: P03231), and gL protein (1-137) (Genebank #: CAD53428).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with gH/gL complex(EBV/B95-8) and gH/gL/gp42 complex (EBV/B95-8)
Storage	Store at 4oC; DO NOT FREEZE; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

ELISA Profile of IT-005-030M9

