



Anti-gp120/gp160(Pan)

Product Name	Anti-gp120/gp160(Pan)
Catalog No	IT-001-002M18
Source	Mouse monoclonal antibody, clone 12C2
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	Recombinant gp120(JRFL)(B) protein (IT-001-0024p)(Genebank No. U63632)
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with most HIV-1 gp120/gp160 proteins as tested. Minimal cross-reactivity to other subtypes.
Storage	Store at 4°C; DO NOT FREEZE ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

ELISA Profile of IT-001-002M18

