

Anti- NP(SFTS Virus)(Huaiyangshan Virus)



Product Name	Anti-NP(SFTS Virus)(Huaiyangshan Virus)
Catalog No	IT-017-006M7
Source	Mouse monoclonal antibody, clone 39C10
Isotype	IgG2a
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing NP (SFTS Virus)(Huaiyangshan Virus) protein (a.a. 1-245) (Genebank Accession No. AEA29891; or ADZ04472).
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts specifically with NP (SFTS Virus)(Huaiyangshan Virus). Minimal cross-reactivity to other Bunyavirus nucleocapsid protein.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

ELISA Profile of IT-017-006M7

