

Anti-Hemagglutinin Protein H (CDV)



Product Name	Anti- Hemagglutinin Protein H (CDV)
Catalog No	IT-007-003M2
Source	Mouse monoclonal antibody, clone 13B2
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing hemagglutinin protein H (Canine distemper virus/strain Onderstepoort)(a.a.144-604)(Gene Bank Accession No.AAK54669)(IT-007-003p)
Applications	Neut., ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with H (Canine distemper virus) protein. Cross-reactivity not determined.
Storage	Stable for 6 months at 4°C from the date of shipment. For long term storage, aliquot and freeze at -80°C.Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

