


HA1(A/environment/California/7451/2010) (H7N3))



Product Name	HA1 (A/environment/California/7451/2010(H7N3))
Catalog No	IT-003-0072p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged HA1 (A/environment/California/7451/2010(H7N3))(a.a.18-349) (See its amino acid sequence below; Genebank #: AFK75166)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; DO NOT FREEZE ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



Protein SEQ:

dkiclgghavangtkvntltergievvnatetvetanikkiciqrkrptdlgqcglgtligppqcdqflefadliierregtdvc
ypgkftneeslrqilrgsggidkepmgftygirtngatsacrrsgssfyamkwllnsdnaafpqmtksyrnprnkpalii
wgvhhsgsateqtklygsgnklitvgsskyqqsfpspgtrpqvngqsggridfhwlldpndvtvftfngafiapdrasflrge
slgiqsdvpldpgcegdchsggtiisslpfqninprtvgkcpvyvkqtslllatgmrnvpenpktrg