

# HA( $\Delta$ TM)(H7N9)(A/Shanghai/2/2013)



---

Product Name	HA( $\Delta$ TM)(H7N9)(A/Shanghai/2/2013)
Catalog No	IT-003-0074 $\Delta$ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 $\mu$ g of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged HA( $\Delta$ TM)(H7N9)(A/Shanghai/2/2013) (a.a.19-524) (GISAID #: EPI439502). This is the HA protein from the outbreak of the novel H7N9 virus in China since late March, 2013.
Applications	Protein standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

