


M1(H1N1) (A/Puerto Rico/8/1934)



Product Name	M1(H1N1)(A/Puerto Rico/8/1934)
Catalog No	IT-003-045Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg of protein in PBS with 8M Urea.
Protein	C-terminal 6xHis tagged M1(H1N1)(A/Puerto Rico/8/1934) matrix protein (a.a. 1-252) (Genebank Accession No. AGK63058).
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>95% , based on SDS PAGE
Storage	Stable for 1 year at -80°C from the date of shipment. Avoiding repeated freeze/thaw. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: M1(H1N1)(A/Puerto Rico/8/1934)