


M1(H7N9) (A/Shanghai/2/2013)



Product Name	M1(H7N9)(A/Shanghai/2/2013)
Catalog No	IT-003-047Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg of protein in PBS with 8M Urea.
Protein	C-terminal 6xHis tagged M1(H7N9)(A/Shanghai/2/2013) matrix protein (a.a. 1-252) (GISAID# EPI439497, same as M1(A/chicken/Taixing/10/2010(H9N2)) that Genbank# is AEO52383).
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>95% , based on SDS PAGE
Storage	Stable for 1 year at -80°C from the date of shipment. Avoiding repeated freeze/thaw. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: M1(H7N9)(A/Shanghai/2/2013)