


Capsid (Chikungunya virus)



Product Name	Capsid (Chikungunya virus)
Catalog No	IT-022-003Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg of protein in PBS with 8M Urea.
Protein	N-terminal GST-tagged and C-terminal 6xHis tagged Chikungunya virus capsid protein (see below for its aa sequence).
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>90% , based on SDS PAGE
Storage	Stable for 6 months at -20°C from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



Capsid Protein Sequence:

efiptqtfnrryqprpwtptrptiqvirprprpqrqagqlaqlisavnkltmravpqqkpararankkqkqkqqapqnntnq
kkqppkkkpaqkkkkpgrrrermcmkiendcifevkhegkvtgyaclvgdkvmkpahvkgtdnadlaklifrssky
dlecaqipvhmksdaskfthekpegyynwhhgavqysggrftiptgagkpgdsgrpifdnkgrvvaivlgganegartal
svvtwnkdivtkitpegae