

Anti-H3(A/Victoria/361/2011)(H3N2)



Product Name	Anti- H3 (A/Victoria/361/2011)(H3N2)
Catalog No	IT-003-00423M10
Source	Mouse monoclonal antibody, clone 31C5
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing HA1(A/Victoria/361/2011)(H3N2) protein (a.a. 18-348)(IT-003-00423p).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Pan H3 as tested. Also reacts with H14. Minimal cross-reactivity to other subtypes.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

