

HA(Δ TM)(B/Massachusetts/2/2012)



Product Name	HA(Δ TM)(B/Massachusetts/2/2012)
Catalog No	IT-003-B19 Δ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged HA(B/Massachusetts/2/2012) (a.a.16-546) (Genbank No. AGL06036)
Applications	Antibody ELISA, Immunogen, WB standard, etc.
Purity	>90%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

