

# Anti-HA1(H7N9)(A/Shanghai/2/2013)



Product Name	Anti-HA1(H7N9)(A/Shanghai/2/2013)
Catalog No	IT-003-0081
Source	Rabbit polyclonal antibody, purified IgG
Formulation	Each vial contains 100 µg of purified IgG (1 mg/ml) in PBS
Immunogen	DNA vaccine expressing HA1 (H7N9)(A/Shanghai/2/2013)(a.a.19-339) (GISAID #: EPI439502). This is the HA1 protein from the outbreak of the novel H7N9 virus in China since late March, 2013.
Applications	WB, ELISA, etc
Specificity	Reacts with HA1(H7N9)(A/Shanghai/2/2013). Cross-reactivity to other HAs not tested.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

## ELISA Profile of IT-003-0081

