

Human IL6 Receptor alpha Chain (Ectodomain)



Product Name	Human IL6 Receptor alpha chain (Ectodomain)
Catalog No	IT-401-002p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged human Human IL6 Receptor alpha chain (Ectodomain) protein (a.a.20-366) (Genebank Accession No. NP_000556)
Applications	Antibody ELISA, Immunogen, WB standard, etc.
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Human IL6 Receptor alpha chain (Ectodomain)