

# Anti-NP (RSV)



Product Name	Anti- NP (RSV)
Catalog No	IT-020-001EM2
Source	Mouse monoclonal antibody, clone 14F8
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	Recombinant NP (RSV)(Human respiratory syncytial virus A2) protein (a.a. 1-391) (Genebank Accession No. P03418). (IT-020-001Ep)
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with NP (RSV)(Human respiratory syncytial virus A2) protein. Cross-reactivity not tested.
Storage	Stable for 6 months at 4°C from the date of shipment. For long term storage, aliquot and freeze at -80°C. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

**ELISA Profile of IT-020-001EM2**

