


# NS1 (DV2)



---

Product Name	NS1 (DV2)
Catalog No	IT-006-0051p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged Dengue Virus-2 NS1 protein + V5 (Genebank No. FJ744720.1).
Applications	WB standard, antibody ELISA, vaccine, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: NS1 (DV2) protein

Note: