


NS1 (WNV)



| | |
|----------------|--|
| Product Name | NS1 (WNV) |
| Catalog No | IT-006-0053p |
| Source | Viral protein purified from 293 cell culture |
| Formulation | Each vial contains 100 µg of purified protein (1 mg/ml) in PBS. |
| Protein | C-terminal 6xHis tagged West Nile virus NS1 protein + V5 (Genebank No. EU081844.1). |
| Applications | WB standard, antibody ELISA, vaccine, etc |
| Purity | >95% |
| Storage | Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |



E: NS1 (WNV) protein

Note: