

# Anti-NS1(Dengue Virus 1)



Product Name	Anti-NS1(DENV1)
Catalog No	IT-006-0050M5
Source	Mouse monoclonal antibody, clone10F5
Isotype	IgG2a
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol
Immunogen	DNA vaccine expressing NS1(Dengue virus 1/GD05/99) (a.a. 1-353) (Genebank No. AAM12552)
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Cross-reacts with NS1 from DENV2, 3, JEV, WNV, TEV.
Storage	Store at 4°C, stable for 6months from the date of shipment. For long term storage, aliquot, and freeze at -80°C. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

## ELISA Profile of IT-006-0050M5

