

Anti-gp120/gp160(Con_A2/2006)



Product Name	Anti- gp120/gp160(Con_A2/2006)
Catalog No	IT-001-001M12
Source	Mouse monoclonal antibody, clone 4B6
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	Recombinant gp120(Con_A2/2006) protein (a.a. 34-518)(IT-001-ConA2p)
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with gp120/gp160(Con_A2/2006) proteins only as tested. Minimal cross-reactivity to other subtypes.
Storage	Store at 4°C; DO NOT FREEZE ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

ELISA Profile of IT-001-001M12

