

# Anti-NS1(Dengue Virus 2)



Product Name	Anti-NS1(DENV2)
Catalog No	IT-006-0051M1
Source	Mouse monoclonal antibody, clone15D1
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol
Immunogen	DNA vaccine expressing NS1(Dengue virus 2/TH-36) (a.a. 1-352) (Genebank No. CAA48937)
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with NS1 from DENV2. Not cross-reactive to NS1 from DENV1, 3, JEV, WNV, and TEV as tested.
Storage	Store at 4°C, stable for 6months from the date of shipment. For long term storage, aliquot, and freeze at -80°C. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

## ELISA Profile of IT-006-0051M1

