

Anti-NS1(Dengue Virus 4)



Product Name	Anti-NS1(DENV4)
Catalog No	IT-006-0055M3
Source	Mouse monoclonal antibody, clone 8C2
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol
Immunogen	DNA vaccine expressing NS1(Dengue virus 4/VE/BID-V2491/2007) (a.a. 1-352) (Genebank No. ACQ44391)
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with NS1 from DENV4 and JEV. Not cross-reactive to NS1 from DENV1, 2, 3, WNV, and TEV as tested.
Storage	Store at 4°C, stable for 6months from the date of shipment. For long term storage, aliquot, and freeze at -80°C. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

ELISA Profile of IT-006-0055M3

