

# Anti-NS1(TEV)



Product Name	Anti-NS1(TEV)
Catalog No	IT-006-0056M3
Source	Mouse monoclonal antibody, clone 10G3
Isotype	IgM
Formulation	Each vial contains 100 µg purified IgM in 100 µl PBS with 40% glycerol
Immunogen	DNA vaccine expressing NS1(Tick-borne encephalitis virus) (a.a. 1-371) (Genebank No. NP_043135)
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with NS1 from TEV. Not cross-reactive to NS1 from DENV1, 2, 3, 4, WNV, JEV and YFV as tested.
Storage	Store at 4°C, stable for 6months from the date of shipment. For long term storage, aliquot, and freeze at -80°C. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

**ELISA Profile of IT-006-0056M3**

