

Anti-NS1(TEV)



| | |
|----------------|---------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name | Anti-NS1(TEV) |
| Catalog No | IT-006-0056M7 |
| Source | Mouse monoclonal antibody, clone 24F12 |
| Isotype | IgG1 |
| Formulation | Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol |
| Immunogen | DNA vaccine expressing NS1(Tick-borne encephalitis virus) (a.a. 1-371) (Genebank No. NP_043135) |
| Applications | WB, ELISA, IP, IF etc. |
| Epitope | Linear |
| Specificity | Reacts with NS1 from TEV, DENV1-3, JEV and WNV.. |
| Storage | Store at 4°C, stable for 6months from the date of shipment. For long term storage, aliquot, and freeze at -80°C. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |

ELISA Profile of IT-006-0056M7

