


VP4(EV71)



Product Name	VP4(EV71)
Catalog No	IT-016-006Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg of lyophilized protein in PBS with 8M Urea. Reconstitute the protein with 100 µl of Millipore water.
Protein	C-terminal GST tagged VP4 region of EV71 (Strain Shanghai 036-2009) capsid protein (a.a.1-69) (GenebankAccession No. ACN59569).
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>95% , based on SDS PAGE
Storage	Before reconstitution, stable for 1 year at -20°C from the date of shipment. After reconstitution, stable for a month at 4°C. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: GST tagged VP4 region of EV71