



Anti-HA(H7N9)(A/Shanghai/2/2013)

Product Name	Anti-HA(H7N9)(A/Shanghai/2/2013)
Catalog No	IT-003-0074M5
Source	Mouse monoclonal antibody, clone 26D8
Isotype	IgG2a
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing HA(H7N9)(A/Shanghai/2/2013) protein (GISAID #: EPI439502)(IT-003-0074p).
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with HA(A/Shanghai/1/2013)(H7N9), HA (A/Shanghai/2/2013)(H7N9), H7(A/Netherlands/03)(H7N7), H7(A/Netherlands/219/03)(H7N7); H8(A/pintail duck/Alberta/114/1979)(H8N4) and H15(A/Australia/83/00)(H15N8) as tested.
Storage	Store at 4°C, stable for 6months from the date of shipment. For long term storage, aliquot, and freeze at -80°C. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

