

# Anti-HA(H7N9)(A/Shanghai/1/2013)



Product Name	Anti- HA(H7N9)(A/Shanghai/1/2013)
Catalog No	IT-003-0073M4
Source	Mouse monoclonal antibody, clone 33B2
Isotype	IgG2a
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing HA(H7N9)(A/Shanghai/1/2013) protein (GISAID #: EPI439486)(IT-003-0073p).
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with HA(A/Shanghai/1/2013)(H7N9), HA (A/Shanghai/2/2013)(H7N9), H7(A/Netherlands/03)(H7N7), H7(A/Netherlands/219/03)(H7N7). No cross-reactivity to other subtype as tested.
Storage	Store at 4°C, stable for 6months from the date of shipment. For long term storage, aliquot, and freeze at -80°C. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

