

Anti-S1(MERS-CoV/EMC 2C/2012)



| | |
|----------------|---|
| Product Name | Anti-S1(MERS-CoV/EMC 2C/2012) |
| Catalog No | IT-002-010M4 |
| Source | Mouse monoclonal antibody, clone 11D12 |
| Isotype | IgG1 |
| Formulation | Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol. |
| Immunogen | DNA vaccine expressing S1 of the spike protein (MERS-CoV/EMC 2C/2012)(Genebank No. AFS88936)(a.a. 18-725)(IT-002-010p). |
| Applications | ELISA, IP, IF etc. |
| Epitope | Conformational |
| Specificity | Reacts with spike protein (MERS-CoV/EMC 2C/2012) only as tested. |
| Storage | Store at 4°C, stable for 6months from the date of shipment. For long term storage, aliquot, and freeze at -80°C. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |

