

# Anti-S1(MERS-CoV/EMC 2C/2012)



Product Name	Anti-S1(MERS-CoV/EMC 2C/2012)
Catalog No	IT-002-010M6
Source	Mouse monoclonal antibody, clone 24B7
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing S1 of the spike protein (MERS-CoV/EMC 2C/2012)(Genebank No. AFS88936)(a.a. 18-725)(IT-002-010p).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with spike protein (MERS-CoV/EMC 2C/2012) only as tested.
Storage	Store at 4°C, stable for 6months from the date of shipment. For long term storage, aliquot, and freeze at -80°C. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

