

Anti-S1(Bat coronavirus HKU4-2)



Product Name	Anti-S1(Bat coronavirus HKU4-2)
Catalog No	IT-002-020M4
Source	Mouse monoclonal antibody, clone 26A9
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing S1 of the spike protein (Bat coronavirus HKU4-2)(Genebank No. ABN10848)(a.a. 21-729)(IT-002-020p).
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with spike protein (Bat coronavirus HKU4-2), and weakly cross-reacts with spike protein (MERS-CoV/EMC 2C/2012).
Storage	Store at 4°C, stable for 6 months from the date of shipment. For long term storage, aliquot, and freeze at -80°C. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

