


Anti-H5(A/Hubei/1/2010)(H5N1)



Product Name	Anti- H5(A/Hubei/1/2010)(H5N1)
Catalog No	IT-003-00519M3
Source	Mouse monoclonal antibody, clone 30A5
Isotype	IgG2a
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing H5 (A/Hubei/1/2010)(H5N1) protein (Genbank No. EPI267032) (IT-003-00519ΔTMp).
Applications	Neut., ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with HA(A/Hubei/1/2010)(H5N1) specifically. No Cross-reactivity to other subtype as tested.
Storage	Store at 4°C, stable for 6 months from the date of shipment. For long term storage, aliquot, and freeze at -80°C. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

ELISA Profile of IT-003-00519M3



ELISA Profile of IT-003-00519M3

