

Antigen MT3413 (M. Tuberculosis/H37Rv)



| | |
|----------------|--|
| Product Name | Antigen MT3413 (M. Tuberculosis/H37Rv) |
| Catalog No | IT-500-008Ep |
| Source | Recombinant protein purified from E.coli |
| Formulation | Each vial contains 100 µg of lyophilized protein in PBS with 8M Urea. Reconstitute the protein with 100 µl of Millipore water. |
| Protein | N-terminal GST fusion protein with C-terminal 6xHis tagged antigen MT3413 (M. Tuberculosis/H37Rv)(a.a.1-107) (Genbank Accession No. NP_337942). |
| Applications | WB standard, antibody ELISA, immunogen, etc. |
| Purity | >95% , based on SDS PAGE |
| Storage | Before reconstitution, stable for 1 year at -20°C from the date of shipment. After reconstitution, stable for a month at 4°C. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |



E: antigen MT3413