

# NS1 (Yellow fever virus)



|                |  |
|----------------|--|
| Product Name   | NS1 (Yellow fever virus)   |
| Catalog No     | IT-006-0057p   |
| Source         | Recombinant protein purified from 293 cell culture   |
| Formulation    | Each vial contains 100 µg of purified protein (1 mg/ml) in PBS   |
| Protein        | C-terminal V5-6xHis tagged NS1 (Yellow fever virus) antigen (a.a.779-1136) (Genbank Accession No. AEQ35299).   |
| Applications   | Antibody ELISA, Immunogen, WB standard, etc.   |
| Purity         | >90% , based on SDS PAGE   |
| Storage        | Before reconstitution, stable for 1 year at -20°C from the date of shipment. After reconstitution, stable for a month at 4°C. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures.   |

