

GP(Sierra-Leone/2014)(ΔTM)



Product Name	GP(Ebola/Sierra-Leone/2014)(ΔTM)
Catalog No	GI-014-005p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged GP(Sierra-Leone/2014)(ΔTM) protein (a.a.33-632) (Genebank #: AIE11800). GP protein from Filoviruses is heavily glycosylated due to its mucin-like-domain.
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

250 kD —
150 kD —
100 kD —
75 kD —
50 kD —
37 kD —
25 kD —
20 kD —
15 kD —
10 kD —

E:GP(Sierra-Leone/2014)(ΔTM)