

# GP(Sierra-Leone/2014) (ΔMUC)(ΔTM)



Product Name	GP(Ebola/Sierra-Leone/2014)(ΔMUC)(ΔTM)
Catalog No	GI-014-006p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis taggedGP(Sierra-Leone/2014)(ΔMUC)(ΔTM)protein (a.a.33-311,463-632) (Genebank #:AIE11800). This protein lacks the mucin-like-domain, but retains the receptor binding domain.
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

250 kD—  
150 kD—  
100 kD—  
75 kD—  
50 kD—  
37 kD—  
25 kD—  
20 kD—  
15 kD—  
10 kD—

E:GP(Sierra-Leone/2014)(ΔMUC)(ΔTM)