

# Anti-GP(Ebola/Sierra-Leone/2014)



---

Product Name	Anti-GP(Ebola/Sierra-Leone/2014)
Catalog No	IT-014-007
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	DNA vaccine expressing ectodomain of glycoprotein of (Ebola/Sierra-Leone/2014).
Applications	ELISA(1:1000-1:100,000), IP, etc.
Specificity	Reacts with GP(Ebola/Sierra-Leone/2014) glycoproteins. Cross-reactivity with other subtypes not determined.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

