

# Anti-NS1(Zika)



|                |  |
|----------------|--|
| Product Name   | Anti-NS1(Zika)   |
| Catalog No     | IT-006-0063M1  |
| Source         | Mouse monoclonal antibody, clone 24F12   |
| Formulation    | Each vial contains 100 µg purified IgG in 100 µl PBS.  |
| Immunogen      | DNA vaccine expressing NS1(JEV) (Genebank No. AJP15451).   |
| Applications   | WB, ELISA(1:1000-1:100,000), IP, etc.  |
| Specificity    | Reacts with NS1(Zika). Cross-reactive with NS1(JEV).   |
| Storage        | Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures.   |

