

Anti-NS1(Zika)



Product Name	Anti-NS1(Zika)
Catalog No	IT-006-0063M3
Source	Mouse monoclonal antibody, clone 51D7
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS.
Immunogen	DNA vaccine expressing NS1(Japanese encephalitis virus) (a.a. 1-382) (Genebank No. YP_006355435).
Applications	WB, ELISA (1:1000-1:100,000), IP, etc.
Specificity	Reacts with NS1(Zika). Cross-reactive with NS1 (JEV).
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

