


NS1 (Zika)



| | |
|----------------|--|
| Product Name | NS1 (Zika) |
| Catalog No | IT-006-0063p |
| Source | Recombinant protein purified from transfected HEK 293 cells |
| Formulation | Each vial contains 100µg protein (1mg/ml) in PBS (pH7.4) |
| Protein | C-terminal 6xHis tagged NS1 (Zika virus, Brazil) antigen (a.a.805- 1147) (Genbank Accession No. ALU33341). |
| Purity | >95% , based on SDS PAGE |
| Storage | Keep refrigerated if use within a month. For long term storage, aliquot it into small volume and store at -80 °C. Avoid repeated freezing and thawing. The product will be expired 6 months after receiving if stored properly. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |



NS1(Zika)