

# Anti-VP1 (EV71)



Product Name	Anti-VP1 (EV71)
Catalog No	IT-016-007M1
Source	Mouse monoclonal antibody, clone 6F5.
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS.
Immunogen	VP1(EV71) antigen (Genebank No. AJP15451).
Applications	WB, ELISA (1:1000-1:10,000), IP, etc.
Specificity	Reacts with VP1 (EV71).
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

