



Anti-gp120(SF162P)(Clade B)

Product Name	Anti- gp120/160(SF162P)(Clade B)
Catalog No	IT-001-145M2
Source	Mouse monoclonal antibody, clone 6E3
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	Recombinant gp140 (SF162P)(Clade B) protein (a.a. 30-668) (Genebank No. AAN05642)
Applications	ELISA, WB, IP, IF etc.
Epitope	Linear
Specificity	Reacts with SHIV gp120/160(SF162P)(Clade B) or related strains only as tested so far. Minimal cross-reactivity to other subtypes.
Storage	Store at 4°C; DO NOT FREEZE ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

