



gp140(SIVsmE660)

Product Name	gp140(SIVsmE660)
Catalog No	IT-001-141p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged gp140(SIVsmE660) protein (a.a.18-700) (Genebank #: ACQ89608)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: gp140(SIVsmE660)

Note: product is partial cleaved at its natural cleavage site