

# Anti-H1N1/H3N2(Seasonal)



Product Name	Anti-H1N1/H3N2(Seasonal)
Catalog No	IT-003-001M11-HRP
Source	Mouse monoclonal antibody, clone 12B9
Isotype	IgG1
Formulation	Each vial contains 50 µg HRP labeled IgG in 100 µl PBS with 40% glycerol
Immunogen	Recombinant HA proteins from seasonal strains
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with most seasonal flu strains. No cross-reactivity to other subtypes.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

