

Anti- HA(H1N1,H3N2(Seasonal))



Product Name	Anti- HA(H1N1,H3N2(Seasonal))
Catalog No	IT-003-001M12
Source	Mouse monoclonal antibody, clone 38D3
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	Sequential immunizations using recombinant H1(ΔTM) (H1N1/A/New Caledonia/20/99)(a.a.18-530)(Genebank No. AAP34324), and H5(H5N1)(A/Vietnam/1203/2004) protein (a.a. 18-530) (GeneBank No. AY651334)
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with seasonal H1 and H3 strains only. Not reactive with 2009 swine flu HA. Minimal cross-reactivity to other subtypes.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

