

# Anti- HA(H1N1(Seasonal))



Product Name	Anti- HA(H1N1(Seasonal))
Catalog No	IT-003-001M13
Source	Mouse monoclonal antibody, clone 85F7
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	Sequential immunizations using recombinant HA1(H1N1) (A/New Caledonia/20/99(H1N1)) ) protein (a.a.18-343) (Genebank No. <a href="#">AAP34324</a> ), and HA1(A/California/06/2009(H1N1)) protein (a.a. 18-344)(Genebank accession <a href="#">ACP41935</a> )
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with seasonal H1 including 2009 swine flu HA. Minimal cross-reactivity to other subtypes.
Storage	Store at 4°C; <b>DO NOT FREEZE</b> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

**ELISA Profile of IT-003-001M13**

