

Anti- HA(H1N1)



Product Name	Anti- HA(H1N1)
Catalog No	IT-003-001M2
Source	Mouse monoclonal antibody
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol
Immunogen	Recombinant HA1 protein (A/New Caledonia/20/99)(Genebank No. ABF21272)
Applications	ELISA, Neut., IP, IF etc.
Epitope	Conformational
Specificity	Reacts with most HA(H1N1) seasonal strains (post 1999). Cross-reactivity to other subtypes not tested.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.