

Anti- HA(H1N1/2009)



Product Name	Anti- HA(H1N1)
Catalog No	IT-003-001M3
Source	Mouse monoclonal antibody
Isotype	IgG1
Formulation	Each vial contains 100µg purified IgG in 100µl PBS with 40% glycerol
Immunogen	Recombinant HA1(A/California/06/2009(H1N1)) protein (a.a. 18-344)(Genebank accession ACP41935).
Applications	ELISA, Neut., IP, IF etc.
Epitope	Conformational
Specificity	Reacts with HA(2009 H1N1) , not react with seasonal strains (post 1957). No cross-reactivity to other subtypes.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

