

Anti- H1(California/06/2009)



Product Name	Anti- H1(California/06/2009(H1N1))
Catalog No	IT-003-001M5-Biotin
Source	Mouse monoclonal antibody, clone 2F3
Isotype	IgG2a
Formulation	Each vial contains 50 µg biotinylated IgG in 200 µl PBS with 40% glycerol
Immunogen	Recombinant HA1(A/California/06/2009(H1N1)) protein (a.a. 18-344)(Genbank accession ACP41935).
Applications	ELISA, Neut., IP, IF etc.
Epitope	Conformational
Specificity	Reacts specifically with H1(2009 H1N1, or swine flu) , not react with seasonal strains. No cross-reactivity to other subtypes.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

