

Anti-H1(A/New Caledonia/20/99)



Product Name	Anti- H1(A/New Caledonia/20/99)
Catalog No	IT-003-001M7
Source	Mouse monoclonal antibody, clone 34-55
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol
Immunogen	Recombinant HA1 protein (A/New Caledonia/20/99)(Genebank No. ABF21272)
Applications	ELISA, Neut., IP, IF, etc.
Epitope	Linear
Specificity	Reacts with H1(A/New Caledonia/20/99). Not cross-react with other H1 and other subtypes.
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

ELISA Profile of IT-003-001M7

