

Anti-H2(H2N2)(A/Canada/720/05)



Product Name	Anti- H2(H2N2)(A/Canada/720/05)
Catalog No	IT-003-003M3
Source	Mouse monoclonal antibody, clone 100E4
Isotype	IgG2a
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol
Immunogen	Recombinant HA ₁ protein (H2N2)(A/Canada/720/05) (Genebank No. AAU28987)
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with most H2(H2N2). Cross-reacts with H6, H5, H8, H13. Not cross-reacts with HA from other subtypes.
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

