

# Anti- H7(H7N7/England/268/96))



Product Name	Anti- H7(A/England/268/96(H7N7))
Catalog No	IT-003-008M7
Source	Mouse monoclonal antibody, clone 17E9
Isotype	IgG2a
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing HA (H7N7/England/268/96) (Genebank No. <a href="#">AF028020</a> ).
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with H7(A/Netherlands/03(H7N7)), H7(A/England/96(H7N7)), H15 (A/Australia/83/00(H15N8)) proteins as tested so far. Minimal cross-reactivity to other subtype.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

**ELISA Profile of IT-003-008M7**

